

# Quality Club™ Specific IgE

Directions for Use 52-5203-EN/32

## INTENDED USE

Quality Club Specific IgE is a quality assessment program for ImmunoCAP Specific IgE users. The individual laboratory's continuous performance is monitored by comparison of test results generated from other ImmunoCAP Specific IgE users.

## PRINCIPLE OF THE PROCEDURE

This kit includes reagents for a period of four months. Each month one of the control samples is assayed and the results are registered online. After a statistical evaluation of all the results, a laboratory specific report is made available.

## REAGENTS AND MATERIAL

The two digit suffix (-xx) on the article number may vary between countries. All kits are not available in all countries.

The expiration date and storage temperature are stated on the labels. Do not use reagents beyond their expiration dates. Make sure to check the expiration date of **all** reagents before use.

### Reagents for Phadia instruments

- **Quality Club Specific IgE** (Art No 10-9298-01: 4 x 0.5 ml)

### Details of reagents

Quality Club Specific IgE	
Control samples Sodium azide 0.05%	Ready for use. Store at 2 – 8 °C until expiration date.
Quality Club Specific IgE Controls are prepared from selected pooled human samples and contain IgE antibodies to certain allergens. The Control Samples contain specific IgE antibodies at varying levels representing three different allergens. Included are major allergens from different groups such as pollens, epithelia, mites and foods.	



### Precautions

- For in vitro diagnostic use. Not for internal or external use in humans or animals.
- The reagents are manufactured from human blood components. The source materials have been tested by immunoassay for hepatitis B surface antigen, for antibodies to HIV1, HIV2 and hepatitis C virus and found to be negative. Nevertheless, all recommended precautions for the handling of blood derivatives should be observed. Please refer to Human Health Service (HHS) Publication No. (CDC) 93-8395 or other local/national guidelines on laboratory safety procedures.

## PROCEDURES

Quality Club control samples should be treated in the same way as patient samples and run in accordance with the Directions for Use valid for the assay used. Make sure to transfer the Quality Club controls to the same type of tubes that are used for the patient samples.

The original Quality Club control vials must not be used in the instrument.

The kit contains four vials of control sample, one vial for each month during the cycle period of four months. The month during which the analysis should be performed is printed on the vial labels. The allergens to be analyzed are printed on each vial.

The results obtained should be registered online before the last day of the month.

### Technical support

If you need any technical support, please contact our local sales office or Phadia AB at qualityclub.idd@thermofisher.com.

### Allergen codes

Allergens used in Quality Club Specific IgE:

d1	Dermatophagoides pteronyssinus
d2	Dermatophagoides farinae
e1	Cat dander
e3	Horse dander
e5	Dog dander
f1	Egg white
f13	Peanut
g3	Dactylis glomerata
g6	Phleum pratense
t3	Betula verrucosa
m3	Aspergillus fumigatus

m6	Alternaria alternata
w6	Artemisia vulgaris

## LIMITATIONS OF THE PROCEDURE

This program is intended for Quality Assessment of the individual laboratory's performance compared to other laboratories and not for control of day-to-day variation.

## WARRANTY

The performance data presented here was obtained using the procedure indicated. Any change or modification in the procedure not recommended by Phadia AB may affect the results, in which event Phadia AB disclaims all warranties expressed, implied or statutory, including the implied warranty of merchantability and fitness for use. Phadia AB and its authorized distributors, in such event, shall not be liable for damages indirect or consequential.

## SYMBOLS

	Use by		In vitro diagnostic medical device
	Batch code		Temperature limitation
	Caution		Consult instructions for use
	Manufacturer		Biological Risks
	Sufficient for		

## Patents/Trademarks

The following designations are trademarks belonging to Phadia AB:

ImmunoCAP, Phadia, Quality Club.

**Trademark change:** Phadia AB has changed the trademarks of the instrument platforms from "UniCAP<sup>®</sup>" and "ImmunoCAP<sup>®</sup>" to "Phadia<sup>®</sup>". The new name has been applied to the instruments and related items, e.g. Software and User Manuals. The trademark "ImmunoCAP<sup>®</sup>" has been removed from the System Reagents. This is a trademark change only; the change has no impact on performance or safety.

## Addresses

**AUSTRIA** Phadia Austria GmbH  
Donau-City-Straße 1  
AT-1220 Vienna  
Tel: +43-1 270 2020 Fax: +43-1 270 202020  
**BELGIUM** Phadia NV/SA  
Pontbeekstraat 2  
BE-1702 GROOT-BIJGAARDEN  
Tel: +32-2 749 55 15 Fax: +32-2 749 55 23  
**BRAZIL** Phadia Diagnósticos Ltda.  
Rua Eugênio de Medeiros, 303 cj 1101C  
05425-000 São Paulo – SP  
Tel: +55 11 2730-3134 Fax: +55 11 2730-3009  
**CHINA** Thermo Fisher Scientific (China) Co., Ltd.  
Building 6-7  
No. 27 Xin Jin Qiao Road  
Shanghai 201102  
P.R. China  
Tel: +86 800 810 5118 Fax: +86 400 650 5118  
**CZECH REPUBLIC** Phadia s.r.o.  
Drahobejlova 1019/27  
19000 PRAHA 9  
Tel: +420 220 518 743 Fax: +420 220 518 743  
**DENMARK** Phadia ApS  
Gydevang 33  
DK-3450 ALLERØD  
Tel: +45-70 23 33 06 Fax: +45-70 23 33 07  
**FINLAND** Phadia Oy  
Ratastie 2  
P.O. Box 100  
FIN-01621 VANTAA  
Tel: +358-9 3291 0110 Fax: +358-9 3291 0531  
**FRANCE** Phadia S.A.S.  
BP 610  
FR-78056 ST QUENTIN-EN-YVELINES CEDEX  
Tel: +33-1 61 37 34 30 Fax: +33-1 30 64 62 37  
**GERMANY** Phadia GmbH  
Postfach 1050 DE-790 10 FREIBURG

Tel: +49-761 47 805-0 Fax: +49-761 47805-338

## HONG KONG

Thermo Fisher Scientific (Hong Kong) Limited  
Unit 11-15, 9/F Tower 1  
Grand Central Plaza  
138 Shatin Rural Committee Road, Shatin  
New Territories, Hong Kong  
Tel: +852 2885 4613 Fax: +852 2567 4447

## INDIA

Phadia India Pvt. Ltd  
Unit No.07, 10 & 11, ground floor,  
Splendor forum, plot no 03  
Jasola, District Centre,  
NEW DEHLI-110025  
Tel: +91 11 461 075 55 / 56 Fax: +91 11 461 075 57

## IRELAND

Phadia Ltd. (Irish Branch)  
27 Oakhill  
Moate  
Co. Westmeath  
Tel: +44 1800 625 167 Fax: +44 1800 625 168  
**ITALY** Phadia S.r.l.  
Via Libero Temolo, 4  
IT-201 26 MILANO  
Tel: +39-0264 163 411 Fax: +39-0264 163 415  
**JAPAN** Phadia K.K.  
NBF Ueno Bldg 9th Floor,  
4-24-11, Higashi -Ueno Taito-ku,  
Tokyo 110-0015, Japan  
Tel: +81-3-5826-1660 Fax: +81-3-5826-1692

## KOREA

Phadia Korea Co. LTD.,  
20 Fl, IT Mirea Tower  
60-21, Gasan-dong Geumcheon-gu  
Seoul 153-801  
Tel: +82-2-2027-5400 Fax: +82-2-2027-5404

## THE NETHERLANDS

Phadia B.V.  
Postbus 696  
NL-3430 AR NIEUWEGEIN  
Tel: +31-30 602 37 00 Fax: +31-30 602 37 09

## NORWAY

Thermo Fisher Scientific  
Phadia AS  
Postboks 4756, Nydalen  
NO-0421 OSLO  
Tel: +47-21 67 32 80 Fax: +47-21 67 32 81

## PORTUGAL

Phadia Sociedade Unipessoal Lda  
Lagoas Park - Edifício n°11 - Piso 0  
PT-2740-270 PORTO SALVO  
Tel: +351-214 23 53 50 Fax: +351-214 21 60 36  
**SOUTH AFRICA** Laboratory Specialities (PTY)  
A Phadia Company  
P.O. Box 1259  
Ferndale 2160

Tel: +27 11 792 6790 Fax: +27 11 793 1064

## SPAIN

Thermo Fisher Scientific  
Phadia Spain SL.  
Avda. Alcalde Barnils nº 70 Planta 2  
Edificio Onada  
08174 Sant Cugat del Vallés, Barcelona  
Tel: +34-935 765 800 Fax: +34-935 765 820

## SWEDEN

Phadia AB,  
Marknadsbolag Sverige  
P O Box 6460 SE-751 37 UPPSALA  
Tel: +46-18 16 60 60 Fax: +46-18 16 63 24

## SWITZERLAND

Phadia AG  
Sennweidstrasse 46  
CH-6312 STEINHAUSEN  
Tel: +41-43 343 40 50 Fax: +41-43 343 40 51  
**TAIWAN** Thermo Fisher Scientific Taiwan Co., Ltd.  
6F-1, No. 85, Jhouzih St., NeiHu District  
Taipei City 11493  
Tel: +886 2 8751 6655 Fax: +886 2 8751 5353

## UNITED KINGDOM

Phadia Ltd  
16 Shenley Pavilions, Chalkdeli Drive  
Shenley Wood, Milton Keynes, MK5 6LB  
Tel: +44-1908 76 91 10 Fax: +44-844 324 94 95  
**USA** Phadia US Inc.  
4169 Commercial Avenue  
Portage, Michigan 49002

Tel: +1 800-346-4364 (Toll Free) Fax: +1 269 492-7541

**OTHER COUNTRIES** Phadia AB,

Distributor Sales

P O Box 6460, SE-751 37 UPPSALA

Tel: +46 18 16 50 00 Fax: +46 18 16 63 65



Phadia AB,  
Rapsgatan 7P, P. O. Box 6460, 751 37 Uppsala, Sweden  
Tel: +46 18 16 50 00 Fax: +46 18 14 03 58



Issued December 2000. Revised April 2016.

© Phadia AB, Uppsala, Sweden.